

Quadriceps tendon tensile test

1. used wet cotton gauze to wrap the ends of sample and glued to both sides of clamp.
2. 1000 precond cycles (except for SR protocol, image acquisition was every 10s. Every 0.1s for stress relaxation)
3. 4.5% max strain.

oks00TR2-QAT-CXXX-06-01 lengths (mm): 9.2285 (for preconditioning),9.408 (fore stress relaxation)

oks00TR2-QAT-CXXX-06-02 lengths (mm): 9.4200 (for preconditioning), 9.6110 (for preconditioning)

oks00TR2-QAT-CXXX-06-03 lengths (mm): 9.609 (for preconditioning), 9.7615 mm (for preconditioning)

