

## Quadreiceps tendon tensile test.

1. used wet cotton gauze to wrap the end of sample and glued to both sides of clamp.
2. 1000 precondition cycles (except for stress relaxation protocol, image acquisition was every 10s. every 0.1s for stress relaxation)
3. 4.5 % max strain
4. 20%/s strain rate

oks00TR2-QAT-CXXX-08-01 lengths (mm): 9.3550 (for preconditioning),9.5755 (for stress relaxation)

oks00TR2-QAT-CXXX-08-02 lengths (mm): 9.5750 (for preconditioning), 9.6145 (for stress relaxation)

oks00TR2-QAT-CXXX-08-03 lengths (mm): 9.689 (for preconditioning), 9.771 (for stress relaxation)

