

## Cartilage 2 day unconfined compression (uc) and confined compression (cc) tests

one uc and one cc test in day 1. cc and uc on day 2 and cc on day 3. 2 hours between tests. sample kept in fridge between tests.


day1: oks00TR6-FMC-MPXX-02-01 (UC): thickness: 1.746 mm: positions: 67.0335, 67.102.


day1: oks00TR6-FMC-MPXX-02-01 (CC): thickness: 1.752 mm mm: positions: 61.3050, 61.3435.

day2: oks00TR6-FMC-MPXX-02-03 (CC): thickness:1.735 mm: positions: 36.9855, 37.0025.

day2: oks00TR6-FMC-MPXX-02-04 (UC): thickness: 1.744 mm: positions: 42.6595, 42.7165.

day 3: oks00TR6-FMC-MPXX-02-05 (CC): thickness: 1.738 mm: positions: 34.0145, 34.0935.





Unconfined compression