

Meniscus unconfined compression test:

1. Preconditioning to 5-15% strain for 1000 cycles
2. Thickness-2.718 mm

oks00TR3-MNS-MPXX-03- 01: find contact->preconditioning->find contact->stress relaxation. Positions: 46.006 (for preconditioning),46.200 (for stress relaxation). Tried to prevent soaking as much as possible before starting test.

3 hours between test

oks00TR3-MNS-MPXX-03-02: recompress to 46.006-> preconditioning-> recompress to 46.200-> stress relax.

